

SHINWA EXPRESS

*Information*2012
November


VOL. 2

ULTRON ES-OVM-C

NEW

The ULTRON ES-OVM-C: Designed to confirm separations of clopidogrel bisulfate as per the USP standards!

Clopidogrel bisulfate (antiplatelet drug)



What's new?

- Highly versatile reversed-phase columns
- Shipping inspections: Clopidogrel bisulfate (USP)
- Minimized lot-to-lot variations and stricter specifications
- Storage solvent: Low viscosity 30% aqueous acetonitrile solution



ULTRON ES-OVM

USP COLUMN LIST L57



The **ULTRON ES-OVM** corresponds to **L57** in the **USP LC column list**.
Can be used for reversed-phase analyses in the same manner as ODS columns.

Citalopram hydrobromide (antidepressant)



HPLC chiral separation applications are available on our website at:
<http://shinwa-cpc.co.jp/en>