

吸着剤

P/N	品名	粒度 (MESH)	容量	Max Temp. (°C)	価格(円)	主な用途
S-130	Shincarbon ST	50~80	10g (約30ml)	330	49,000	無機ガス、 C ₁ ~C ₂ 有機ガス
A-1	Activated Alumina	30~60	100ml	250	9,000	有機ガス、水素ガス
A-2		60~80			9,000	
A-2D		80~100			9,000	
A-3	Activated Charcoal	30~60	100ml	220	8,800	有機ガス、水素ガス
A-4		60~80			8,800	
A-4D		80~100			8,800	
M-1	Molecular Sieve 5A	30~60	100ml	250	9,200	有機ガス、メタンなど
M-2		60~80			10,200	
M-2D		80~100			10,200	
M-3	Molecular Sieve 13X	30~60	100ml	250	9,000	有機ガス、メタンなど
M-4		60~80			10,200	
M-4D		80~100			10,200	
S-1	Silica Gel	40~70	100ml	220	8,800	有機ガス、有機ガスなど
S-2		60~80			8,800	
S-2D		80~100			8,800	
	Molecular Sieve 5A	1/8 inch	500g	250	11,100	
		1/16 inch			11,100	
	Molecular Sieve 13X	1/8 inch	500g	250	11,800	
		1/16 inch			11,800	
	Molecular Sieve 3A	1/8 inch	500g	250	11,100	
		1/16 inch			11,100	
	Molecular Sieve 4A	1/8 inch	500g	250	11,100	
		1/16 inch			11,100	